



Product Datasheet

ABRA1 antibody (orb44678)

Description	Mouse Monoclonal to ABRA1.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IP, WB
Immunogen	Recombinant protein corresponding to amino acids 1-313 of ABRA1 with N-terminal His6 tag
Target	ABRA1
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-2 µg/ml. Intracellular staining. Western blotting: Reducing conditions.
Isotype	Mouse IgG2a kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q6UWZ7
Entrez	84142
Dilution Range	Flow cytometry: Recommended dilution: 1-2 µg/ml. Intracellular staining. Western blotting: Reducing conditions.
Expiration Date	12 months from date of receipt.



Western blotting detection of ABRA1-EGFP...



Separation of MOLT-4 cells stained using...