


Product Datasheet

CD26 antibody (orb44682)

Description	Mouse monoclonal antibody to CD26
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	A human T cell clone
Target	CD26
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-2 µg/ml; positive control: human peripheral blood leukocytes.
Isotype	Mouse IgG2a
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P27487
Entrez	1803
Dilution Range	Flow cytometry: Recommended dilution: 1-2 µg/ml; positive control: human peripheral blood leukocytes.
Expiration Date	12 months from date of receipt.



Flow cytometry (surface staining) of hum...