



Product Datasheet

alpha 2 Macroglobulin antibody (orb44694)


Description	alpha 2 Macroglobulin antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Alpha-2-macroglobulin protein (1216-1467AA)
Target	A2M
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P01023
Dilution Range	IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human en...



Immunofluorescence analysis of HepG2 cel...



Immunohistochemical staining of human en...