



Product Datasheet

CD107a antibody (orb44715)

Description	Mouse Monoclonal to CD107a.
Reactivity	Human, Mouse, Primate
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-Fr, IHC-P, WB
Immunogen	Human PBMC
Target	CD107a
Preservatives	Phosphate buffered saline (PBS), pH 7.4
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Western blotting: Recommended dilution: 1 µg/ml. Immunohistochemistry: Recommended dilution: 2-8 µg/ml. Flow cytometry: Intracellular and extracellular staining; recommended dilution: 1-12 µg/ml.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P11279
Entrez	3916
Dilution Range	Western blotting: Recommended dilution: 1 µg/ml. Immunohistochemistry: Recommended dilution: 2-8 µg/ml. Flow cytometry: Intracellular and extracellular staining; recommended dilution: 1-12 µg/ml.
Expiration Date	12 months from date of receipt.



Flow cytometry analysis (intracellular staining) of CD107a in PBMCs.



Immunohistochemistry staining of CD107a in tissue sections.