



Product Datasheet

CD28 antibody (orb44737)

Description	Mouse Monoclonal to CD28.
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	FC, IP
Immunogen	Mouse A20J B lymphoma cells transfected with rat CD28
Target	CD28
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 3-12 µg/ml.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P31042
Entrez	25660
Dilution Range	Flow cytometry: Recommended dilution: 3-12 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of rat CD28 positive splenocytes...



Flow cytometry surface staining pattern ...