



Product Datasheet

CD361 antibody (orb44740)

Description	Mouse Monoclonal to CD361.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IP
Immunogen	Raji cells
Target	CD361
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml, Positive control, Positive control: Raji, Daudi, HL-60 cells, peripheral blood lymphocytes (strongly positive on CD19+ cells), negative control: Jurkat, U-937 cells.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P34910
Entrez	2124
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml, Positive control, Positive control: Raji, Daudi, HL-60 cells, peripheral blood lymphocytes (strongly positive on CD19+ cells), negative control: Jurkat, U-937 cells.
Expiration Date	12 months from date of receipt.



Separation of human neutrophil granulocy...



Flow cytometry surface staining pattern ...