

Product Datasheet

IKK alpha/beta (phospho-Ser176/Ser180) antibody (Biotin) (orb447436)

Description	IKK alpha/beta (phospho-Ser176/Ser180) antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Biotin
Tested Applications	ELISA, ICC, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human IKK alpha/beta around the phosphorylation site of Ser176/Ser180
Target	Phospho-IKK alpha/beta (Ser176 + Ser180)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O15111
Dilution Range	IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.