


Product Datasheet

NFATC1 antibody (orb44745)

Description	NFATC1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human Nuclear factor of activated T-cells, cytoplasmic 1 protein (48-406AA)
Target	NFATC1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O95644
Dilution Range	IHC-P:1:1000-2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human pa...