



Product Datasheet

GFP antibody (orb44764)

Description	Rabbit Polyclonal to GFP.
Reactivity	All
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	EGFP, a native full-length protein
Target	GFP
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Western blotting: Recommended dilution 0.5-1.5 µg/ml; positive control: transfected cells; negative control: non-transfected cells. Immunocytochemistry: Recommended dilution 1-3 µg/ml.
Isotype	Rabbit polyclonal
Clonality	Monoclonal
Purity	Purified by ligand affinity chromatography.
Uniprot ID	P42212
Dilution Range	Western blotting: Recommended dilution 0.5-1.5 µg/ml; positive control: transfected cells; negative control: non-transfected cells. Immunocytochemistry: Recommended dilution 1-3 µg/ml.
Expiration Date	12 months from date of receipt.



Immunoprecipitation of GFP-NLS from HEK293T cells



Immunocytochemistry staining (confocal microscopy) of cells transfected with EGFP