


Product Datasheet

Kinesin antibody (orb44767)

Description	Rabbit Polyclonal to Kinesin.
Reactivity	Human, Mouse, Porcine
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	stalk domain of human kinesin (aa 331-906) expressed in E. coli (FKHC3)
Target	Kinesin
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Isotype	Rabbit polyclonal
Clonality	Monoclonal
Purity	Purified by ligand affinity chromatography.
Expiration Date	12 months from date of receipt.



Western blotting detection of kinesin in...