


Product Datasheet

SLP76 antibody (orb44776)

Description	Mouse Monoclonal to SLP76.
Reactivity	Human, Mouse, Porcine
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Bacterially expressed fusion protein representing amino acids 216-434 of human SLP76 with histidine tag
Target	SLP76
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Immunohistochemistry (paraffin sections): Recommended dilution: 10 µg/ml; positive tissue: thymus
Isotype	Mouse IgG2b
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q13094
Entrez	3937
Dilution Range	Immunohistochemistry (paraffin sections): Recommended dilution: 10 µg/ml; positive tissue: thymus
Expiration Date	12 months from date of receipt.



Western blotting analysis (reducing cond...