

Product Datasheet

SLP76 antibody (orb44777)

Description

Rabbit Polyclonal to SLP76.


Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

Bacterially expressed fusion protein representing amino acids 216-434 of human SLP76 with histidine tag

Immunohistochemistry staining of human t...

Target

SLP76

Preservatives

Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Concentration

1 mg/ml

Storage

Store at 2-8°C. Do not freeze.

Note

For research use only

Application notes

Western blotting: Positive control: JURKAT human T cell leukemia cell lysate, reducing conditions;
Immunohistochemistry (paraffin sections):
Recommended dilution: 10 µg/ml, positive tissue: thymus.

Isotype

Rabbit polyclonal

Clonality

Monoclonal

Purity

Purified by protein-A affinity chromatography.

Uniprot ID
Q13094
Entrez
3937
Dilution Range

Western blotting: Positive control: JURKAT human T cell leukemia cell lysate, reducing conditions;
Immunohistochemistry (paraffin sections):
Recommended dilution: 10 µg/ml, positive tissue: thymus.

Expiration Date

12 months from date of receipt.