

Product Datasheet

atpB antibody (Biotin) (orb44802)

Description	atpB antibody (Biotin)
Species/Host	Rabbit
Reactivity	A. thaliana
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Arabidopsis thaliana ATP synthase subunit beta, chloroplastic protein (92-393AA)
Target	atpB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P19366
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.