

Product Datasheet

SRGAP2C antibody (Biotin) (orb44814)

Description	SRGAP2C antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human SLIT-ROBO Rho GTPase-activating protein 2C protein (1-242AA)
Target	SRGAP2C
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	PODJJO
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.