



Product Datasheet

CD86 antibody (orb44816)

Description	Rat Monoclonal to CD86.
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FA, FC, ICC, IHC-Fr, IP
Immunogen	LPS-activated CBA/Cs mouse splenic B cells
Target	CD86
Preservatives	Phosphate buffered saline (PBS), pH 7.4
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Low Endotoxin, Azide free format. Product does not contain preservatives and is intended primarily for functional applications. Preservative containing buffers are available for IHC applications
Isotype	Rat IgG2a kappa
Clonality	Monoclonal
Purity	Purified by protein-G affinity chromatography.
Uniprot ID	P42082
Entrez	12524
Dilution Range	Functional application: Blocking. Flow cytometry: Recommended dilution: 2 µg/ml; positive control: murine splenocytes.
Expiration Date	12 months from date of receipt.



Separation of murine CD86 positive myelo...



Flow cytometry surface staining pattern ...