



Product Datasheet

CD86 antibody (FITC) (orb44818)

Description	Rat Monoclonal to CD86.
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	LPS-activated CBA/Cs mouse splenic B cells
Target	CD86
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	0.5 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-2 µg/ml; positive control: murine splenocytes.
Isotype	Rat IgG2a kappa
Clonality	Monoclonal
Purity	Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	P42082
Entrez	12524
Dilution Range	Flow cytometry: Recommended dilution: 1-2 µg/ml; positive control: murine splenocytes.
Expiration Date	12 months from date of receipt.



Flow
cytometry
surface
staining
pattern ...



Separation of
murine CD86
positive
myelo...