



Product Datasheet

CD11b antibody (orb44821)

Description	Rat monoclonal antibody to CD11b
Reactivity	Human, Mouse, Primate, Rabbit
Conjugation	Unconjugated
Tested Applications	FA, FC, IHC-Fr, IP
Immunogen	B10 mouse spleen cells enriched for T cells
Target	CD11b
Preservatives	Phosphate buffered saline (PBS), pH 7.4
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Low Endotoxin, Azide free format. Product does not contain preservatives and is intended primarily for functional applications. Preservative containing buffers are available for IHC applications
Isotype	Rat IgG2b
Clonality	Monoclonal
Purity	Purified by protein-G affinity chromatography.
Uniprot ID	G5E8F1
Entrez	16409
Dilution Range	Functional application: In vitro blocking of CD11b. Flow cytometry: Recommended dilution: 1 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of murine CD11b positive myel...



Flow cytometry surface staining pattern ...