



Product Datasheet

CD11b antibody (FITC) (orb44822)

Description	Rat monoclonal antibody against CD11b conjugated to FITC
Reactivity	Human, Mouse, Primate, Rabbit
Conjugation	FITC
Tested Applications	FC
Immunogen	B10 mouse spleen cells enriched for T cells
Target	CD11b
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	0.5 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-3 µg/ml.
Isotype	Rat IgG2b
Clonality	Monoclonal
Purity	Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	G5E8F1
Entrez	16409
Dilution Range	Flow cytometry: Recommended dilution: 1-3 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of
mouse
CD11b
positive
myelo...



Flow
cytometry
surface
staining
pattern ...