


Product Datasheet

CD11a antibody (FITC) (orb44826)

Description	Rat Monoclonal to CD11a.
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	C57BL/6 mouse splenic secondary cytotoxic T lymphocytes
Target	CD11a
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	0.5 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-2 µg/ml.
Isotype	Rat IgG2a
Clonality	Monoclonal
Purity	Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	E9Q5M7
Entrez	16408
Dilution Range	Flow cytometry: Recommended dilution: 1-2 µg/ml.
Expiration Date	12 months from date of receipt.



Flow cytometry analysis (surface stainin...