

Product Datasheet

Neurotensin Receptor 1 antibody (Biotin) (orb448464)

Description	Neurotensin Receptor 1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human NTR1/Neurotensin Receptor 1
Target	Neurotensin Receptor 1
Preservatives	Liquid in PBS containing 50% glycerol, 0.5% rAlbumin and 0.32% sodium azide.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P30989
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.