

Product Datasheet

**GluA4/Ionotropic Glutamate receptor 4 antibody
(Biotin) [Out of stock] (orb448470)**

Description

Rabbit polyclonal antibody against

receptor 4 conjugated to Biotin

Species/Host

Rabbit

Conjugation

Biotin

Immunogen

KLH conjugated synthetic peptide derived from human GluA4/Ionotropic Glutamate receptor 4

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

98kDa

Uniprot ID
P19493, P48058, Q9Z2W8
Dilution Range

WB=1:50-200 ELISA=1:100-1000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.

Expiration Date

12 months from date of receipt.