

Product Datasheet

MGLUR1A/Metabotropic Glutamate Receptor 1a antibody (Biotin) (orb448471)

Description	MGLUR1A/Metabotropic Glutamate Receptor 1a antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human, Mouse, Zebrafish
Conjugation	Biotin
Tested Applications	ELISA, ICC, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human MGLUR1
Target	GRM1
Preservatives	0.01M PBS(pH7.4)
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13255
Dilution Range	IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.