

Product Datasheet

Mouse Anti-Human IgG (Fc) antibody (HRP) (orb44868)

Description

Mouse Monoclonal to Human IgG (Fc).

Reactivity

Human

Conjugation

HRP

Tested Applications

ELISA, WB

Immunogen

Fusion protein of human IgG Fc fragment.

Target

IgG (Fc)

Preservatives

Stabilizing phosphate buffered saline (PBS), pH 7.4, 0.01% thimerosal

Concentration

1 mg/ml

Storage

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Note

For research use only

Application notes

Western blotting and ELISA: The peroxidase conjugate of this antibody is suitable for detection of human IgG Fc fragments. The antibody EM-07 does not crossreact with IgM.

Isotype

Mouse IgG1

Clonality

Monoclonal

Purity


Purified antibody is conjugated with activated horseradish peroxidase (HRP) of high specific enzyme activity under optimum conditions and unconjugated antibody and free HRP are removed by size-exclusion chromatography.

Dilution Range

Western blotting and ELISA: The peroxidase conjugate of this antibody is suitable for detection of human IgG Fc fragments. The antibody EM-07 does not crossreact with IgM.

Expiration Date

12 months from date of receipt.



ELISA
detection of
IgG in
human
serum us...