


Product Datasheet

Mouse Anti-Human IgM antibody (orb44872)

Description	Mouse Monoclonal to Human IgM.
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	IgM from human serum.
Target	IgM
Preservatives	Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Isotype	Mouse IgG2b
Clonality	Monoclonal
Purity	Purified by sequential steps of physicochemical fractionation (differential precipitation and solid-phase chromatography methods).
Expiration Date	12 months from date of receipt.



ELISA analysis of IgM in human serum using...