

Product Datasheet

CSF3R antibody (orb44878)

Description	Goat polyclonal antibody to CSF3R
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Target	G-CSF receptor
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Concentration	0.5 mg/ml
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Preliminary experiments in Human Placenta, Bone marrow and peripheral lymphocytes lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results bee
Protein Sequence	TQGDSILDCVPKDGQ
Clonality	Polyclonal
Entrez	1441
Dilution Range	ELISA: 1:4000, WB: 1 µg/ml
Expiration Date	12 months from date of receipt.