


Product Datasheet

NRXN1 antibody (orb44881)

Description	NRXN1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Neurexin-1-beta protein (184-363AA)
Target	NRXN1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P58400
Dilution Range	WB:1:2000-5000, IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HepG2 cel...