

Product Datasheet

C8ORF34 antibody (Biotin) (orb449260)

Description	C8ORF34 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human C8ORF34
Target	C8ORF34
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q49A92
Expiration Date	12 months from date of receipt.