

Product Datasheet

HIF Prolyl Hydroxylases antibody (Biotin) (orb449379)

Description	HIF Prolyl Hydroxylases antibody (Biotin)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA, ICC, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human HIF Prolyl Hydroxylases
Target	HIF Prolyl Hydroxylases
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NXG6
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.