

Product Datasheet

DCAF8 antibody (Biotin) (orb45040)

Description	DCAF8 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human DDB1- and CUL4-associated factor 8 protein (1-135AA)
Target	DCAF8
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q5TAQ9
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.