

Product Datasheet

RABGAP1L antibody (Biotin) (orb45060)

Description	RABGAP1L antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Rab GTPase-activating protein 1-like, isoform 10 protein (11-262AA)
Target	RABGAP1L
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	B7ZAP0
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.