

Product Datasheet

**Kinesin 5C antibody (Biotin) [Out of stock]
(orb450664)**

Description

Rabbit polyclonal antibody against

Species/Host

Rabbit

Conjugation

Biotin

Tested

ELISA, ICC, IF, IHC-Fr, IHC-P, WB

Applications
Immunogen

KLH conjugated synthetic peptide
derived from human Kinesin
5C/NKHC2

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin,
0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid
repeated freeze/thaw cycles. The
lyophilized antibody is stable at room
temperature for at least one month
and for greater than a year when kept
at -20°C. When reconstituted in sterile
pH 7.4 0.01M PBS or diluent of
antibody the antibody is stable for at
least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

109kDa

Uniprot ID
O60282, **P56536**
Dilution Range

WB: 1:50-200, ELISA: 1:100-1000,
IHC-P: 1:50-200, ICC: 1:50-200, IF:
1:50-200

Expiration Date

12 months from date of receipt.