

Product Datasheet

hydroxyproline antibody (Biotin) (orb451098)

Description

Rabbit polyclonal antibody against


Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

Biotin

Tested

ELISA, IHC-Fr, IHC-P

Applications

Immunohistochemical staining of rat brai...

Immunogen

KLH conjugated synthetic peptide derived from human hydroxyproline

Target

51-35-4

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid : 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.


Immunohistochemical staining of human lu...

Concentration

1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.



Immunohistochemical staining of rat brai...

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Affinity purified by Protein A

CAS Number

51-35-4

Dilution Range

IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000

Expiration Date

12 months from date of receipt.