

Product Datasheet

DTX4 antibody (HRP) (orb45146)

Description	DTX4 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human E3 ubiquitin-protein ligase DTX4 protein (167-432AA)
Target	DTX4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y2E6
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.