

Product Datasheet

Syntaxin 1A + Syntaxin 1B antibody (Biotin) [Out of stock] (orb452649)

Description

Rabbit polyclonal antibody against

conjugated to Biotin

Species/Host

Rabbit

Conjugation

Biotin

Tested

ELISA, ICC, IF, IHC-Fr, IHC-P, WB

Applications
Immunogen

KLH conjugated synthetic peptide derived from human Syntaxin 1A + Syntaxin 1B

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

- 100 µg (in 100 µl): 1 mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

33kDa

Uniprot ID
P61265, P61264, P61268, P61266, P61267
Dilution Range

WB: 1:50-200, ELISA: 1:100-1000, IHC-P: 1:50-200, ICC: 1:50-200, IF: 1:50-200

Expiration Date

12 months from date of receipt.