

## Product Datasheet

### Tubulin Alpha antibody (Biotin) (orb452727)

|                            |                                                                                  |
|----------------------------|----------------------------------------------------------------------------------|
| <b>Description</b>         | Tubulin Alpha antibody (Biotin)                                                  |
| <b>Species/Host</b>        | Rabbit                                                                           |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                |
| <b>Conjugation</b>         | Biotin                                                                           |
| <b>Tested Applications</b> | ELISA, ICC, IHC-Fr, IHC-P, WB                                                    |
| <b>Immunogen</b>           | KLH conjugated synthetic peptide derived from human TUBA1A + TUBA1B              |
| <b>Target</b>              | Tubulin alpha (Loading Control)                                                  |
| <b>Preservatives</b>       | Preservative: 15mM Sodium Azide,<br>Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4 |
| <b>Form/Appearance</b>     | Liquid                                                                           |
| <b>Concentration</b>       | 1mg/ml                                                                           |
| <b>Storage</b>             | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                            |
| <b>Isotype</b>             | IgG                                                                              |
| <b>Clonality</b>           | Polyclonal                                                                       |
| <b>Uniprot ID</b>          | <a href="#">Q71U36</a>                                                           |
| <b>Dilution Range</b>      | WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100,ELISA=1:5000-10000       |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                  |