

Product Datasheet

Bcl-2 antibody (Biotin) (orb452752)

Description	Bcl-2 antibody (Biotin)
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	Recombinant human Bcl-2 protein
Target	Bcl-2
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	26kDa
Uniprot ID	P49950
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.