

Product Datasheet

P4HB antibody (Biotin) (orb45289)

Description	P4HB antibody (Biotin)
Species/Host	Rabbit
Reactivity	Hamster
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Cricetulus griseus Protein disulfide-isomerase protein (20-509AA)
Target	P4HB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8R4U2
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.