

Product Datasheet

GAPDH-Loading Control antibody (Biotin) (orb452926)

Description	GAPDH-Loading Control antibody (Biotin)
Species/Host	Mouse
Reactivity	Hamster, Human, Mouse, Porcine, Rat
Conjugation	Biotin
Tested Applications	ICC, IHC-P, WB
Immunogen	Recombined Human GAPDH
Target	GAPDH-Loading Control
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P04406
Dilution Range	WB=1:10000-1:100000,IHC-P=1:100- 500,ICC=1:100
Expiration Date	12 months from date of receipt.