

Product Datasheet

**alpha Synuclein (phospho-Ser129) antibody
(Biotin) [Out of stock] (orb453138)**

Description

Rabbit polyclonal antibody

(phospho-Ser129) conjugated to Biotin

Species/Host Rabbit

Conjugation Biotin

Tested Applications ELISA, IHC-P, WB

Immunogen Please contact us for immunogen information

Target alpha Synuclein (phospho-Ser129) (Biotin)

Preservatives 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance Liquid: 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Concentration - 100 µg (in 100 µl): 1 mg/ml

Storage Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note For research use only

Isotype IgG

Clonality Polyclonal

Purity Affinity purified by Protein A

MW 15 kDa

Uniprot ID [Q16143](#)

Dilution Range WB: 1:50-200, ELISA: 1:100-1000, IHC-P: 1:50-200

Expiration Date 12 months from date of receipt.