

Product Datasheet

Acetyl NF KB P65 (Lys314 antibody (Biotin) (orb453327)

Description	Acetyl NF KB P65 (Lys314 antibody (Biotin)
Species/Host	Mouse
Reactivity	Mouse, Rat
Conjugation	Biotin
Tested Applications	ICC, IHC-Fr, IHC-P
Immunogen	KLH conjugated Synthesised acetylpeptide derived from human NFKBp65 around the acetylation site of Lys314/315
Target	NF kB P65 (Acetyl-Lys314/Lys315)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	Q04206
Dilution Range	IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.