

Product Datasheet

AVL9 antibody (Biotin) (orb45353)

Description	AVL9 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Late secretory pathway protein AVL9 homolog protein (162-432AA)
Target	AVL9
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NBF6
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.