

Product Datasheet

Semaphorin 5A antibody (Biotin) (orb453562)

Description	Semaphorin 5A antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human Semaphorin 5A
Target	Semaphorin 5A
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13591
Expiration Date	12 months from date of receipt.