

Product Datasheet

GSK-3 beta antibody (Biotin) (orb453620)

Description	GSK-3 beta antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human GSK-3 Beta
Target	GSK-3 beta
Preservatives	0.01M PBS, 0.5%rAlbumin, 0.03% Proclin300
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P49841
Dilution Range	WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.