

Product Datasheet

LC3A/B + LC3B + LC3C antibody (Biotin) (orb453766)

Description	LC3A/B + LC3B + LC3C antibody (Biotin)
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Biotin
Tested Applications	ELISA, ICC, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human LC3A/B + LC3B + LC3C
Target	LC3A/B + LC3B + LC3C
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BXW4
Dilution Range	IHC-P=1:100-500,IHC-F=1:100- 500,ICC=1:100-500,ELISA=1:5000- 10000
Expiration Date	12 months from date of receipt.