

Product Datasheet

GALNT11 antibody (HRP) (orb45383)

Description	GALNT11 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Polypeptide N-acetylgalactosaminyltransferase 11 protein (420-606AA)
Target	GALNT11
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NCW6
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.