

Product Datasheet

IL4I1 antibody (FITC) (orb45392)

Description	IL4I1 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human L-amino-acid oxidase protein (139-339AA)
Target	IL4I1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96RQ9
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.