



Product Datasheet

EXD2 antibody (orb45518)

Description	EXD2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Exonuclease 3\\'-5\\' domain-containing protein 2 protein (137-245AA)
Target	EXD2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NVH0
Dilution Range	WB:1:2000-10000, IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human gastric tissue



Immunohistochemical staining of human testis tissue