

Product Datasheet

PI3 Kinase beta antibody (FITC) [Out of stock] (orb4552)

Description	Rabbit polyclonal to PI3 Kinase beta. Phosphorylates PtdIns, PtdIns4P and PtdIns(4,5)P2 with a prefe...
Species/Host	Rabbit
Conjugation	FITC
Immunogen	KLH conjugated synthetic peptide derived from human PI3 Kinase p110 beta (701-800aa)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	- 100 µg (in 100 µl): 1 mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	123kDa
Uniprot ID	Q9Z1L0 , Q8BTI9 , P42338 , Q3U4Q1 , Q24JU2
Expiration Date	12 months from date of receipt.