

Product Datasheet

TDRD1 antibody (Biotin) (orb45533)

Description	TDRD1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Tudor domain-containing protein 1 protein (12-314AA)
Target	TDRD1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BXT4
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.