



Product Datasheet

SMOC1 antibody (orb45538)


Description	SMOC1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human SPARC-related modular calcium-binding protein 1 protein (277-382AA)
Target	SMOC1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid: Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	> 95%, Protein G purified
Uniprot ID	Q9H4F8
Dilution Range	WB:1:1000-5000, IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse skeletal ...



Immunohistochemical staining of human CO...



Immunofluorescence analysis of HepG2 cel...